

For Research Use Only

PRKDC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61098-3-PBS



Basic Information

Catalog Number:

61098-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29173

GenBank Accession Number:

NM_006904

GeneID (NCBI):

5591

UNIPROT ID:

P78527

Full Name:

protein kinase, DNA-activated,
catalytic polypeptide

Calculated MW:

469 kDa

Purification Method:

Protein G purification

CloneNo.:

3F6E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51771-1, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-3-PBS.



Cytometric bead array standard curve of MP51771-1, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-3-PBS. Standard:Ag29173. Range: 0.156-20 ng/mL.