

For Research Use Only

MAP3K11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61102-1-PBS



Basic Information

Catalog Number: 61102-1-PBS	GenBank Accession Number: BC011263	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 4296	CloneNo.: 1D4C8
Isotype: IgG1	UNIPROT ID: Q16584	
Immunogen Catalog Number: AG29748	Full Name: mitogen-activated protein kinase kinase kinase 11	
	Calculated MW: 847 aa, 93 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51737-1, MAP3K11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61102-1-PBS. Detection antibody: 61102-2-PBS.



Cytometric bead array standard curve of MP51737-1, MAP3K11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61102-1-PBS. Detection antibody: 61102-2-PBS. Standard: Ag29748. Range: 0.156-20 ng/mL