

For Research Use Only

MAP3K11 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61102-2-PBS



Basic Information

Catalog Number:

61102-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29748

GenBank Accession Number:

BC011263

GeneID (NCBI):

4296

UNIPROT ID:

Q16584

Full Name:

mitogen-activated protein kinase
kinase kinase 11

Calculated MW:

847 aa, 93 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G1C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51737-1, MAP3K11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61102-1-PBS. Detection antibody: 61102-2-PBS.



Cytometric bead array standard curve of MP51737-1, MAP3K11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61102-1-PBS. Detection antibody: 61102-2-PBS. Standard: Ag29748. Range: 0.156-20 ng/mL