

For Research Use Only

# CXCR3 Monoclonal antibody

Catalog Number: 65008-1-Ig



## Basic Information

Catalog Number:

65008-1-Ig

Size:

1330 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

AF469635

GeneID (NCBI):

2833

UNIPROT ID:

P49682

Full Name:

chemokine (C-X-C motif) receptor 3

Calculated MW:

46 kDa

Observed MW:

35-40 kDa, 38 kDa, 43 kDa

Purification Method:

Protein G purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data