

For Research Use Only

CPS1 Monoclonal antibody

Catalog Number: 65011-1-Ig **2 Publications**



Basic Information

Catalog Number:

65011-1-Ig

Size:

670 µg/ml

Source:

Mouse

Isotype:

IgG

GenBank Accession Number:

BC020695

GeneID (NCBI):

1373

UNIPROT ID:

P31327

Full Name:

carbamoyl-phosphate synthetase 1,
mitochondrial

Calculated MW:

165 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : Human liver,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jana Deitersen	34059630	Cell Death Dis	WB
Ke Wang	39011520	Biomater Res	WB

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

