

For Research Use Only

# Anti-Mouse CSF1R/CD115 (AFS98)

Catalog Number: 65095-1-Ig



## Basic Information

Catalog Number:

65095-1-Ig

Size:

500ug, 500 µg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug Anti-Mouse CD115 (c-fms) (AFS98) (65095-1-Ig, Clone: AFS98) or 0.5 ug Rat IgG2a Isotype Control (2A3) (65209-1-Ig, Clone: 2A3) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). Cells were not fixed.