

Anti-Mouse IL-12/IL-23 p40 (C17.8)

Catalog Number: 65097-1-Ig

Basic Information

Catalog Number: 65097-1-Ig	GenBank Accession Number: BC103608	Purification Method: Affinity purification
Size: 100ug, 500 µg/ml	GeneID (NCBI): 16160	CloneNo.: C17.8
Source: Rat	UNIPROT ID: P43432	
Isotype: IgG2a, kappa	Full Name: interleukin 12b	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
Mouse

Background Information


Interleukin (IL)-12 and IL-23 are heterodimeric cytokines that share a common p40 subunit (PMID: 11114383). IL-12 is composed of the IL-12 p40 subunit linked to the IL-12 p35 subunit, and the heterodimer signals through the IL-12 receptor (IL-12R), which comprises the IL-12R β 1 and IL-12R β 2 subunits. IL-23 is composed of the IL-23 p19 subunit and the IL-12 p40 (IL-12/23p40) subunit, which signals through IL-23R and IL-12R β 1 (PMID: 11114383; 26121196). IL-12/IL-23 p40 also exists as a monomer and as a homodimer which can act as a potent IL-12 antagonist (PMID: 8958912; 18783467). IL-12/IL-23 p40 is produced by antigen-presenting cells, such as dendritic cells (DCs), monocytes, macrophages, neutrophils and, to a lesser extent, B cells (PMID: 20476918).

Storage

Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

Selected Validation Data



1X10⁶ untreated or IFN γ , LPS, BFA and Monensin treated RAW 264.7 cells were intracellularly stained with 0.5 ug Anti-Mouse IL-12/IL-23 p40 (65097-1-Ig, Clone: C17.8) and FITC anti-rat IgG2a Antibody at dilution 1:100. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).