

For Research Use Only

Anti-Mouse CD275 (B7-H2) (HK5.3)

Catalog Number: 65114-1-Ig



Basic Information

Catalog Number:

65114-1-Ig

Size:

100ug,

Source:

Rat

Isotype:

Rat IgG2a, kappa

GenBank Accession Number:

BC029227

GeneID (NCBI):

50723

Full Name:

icos ligand

Purification Method:

Affinity purification

CloneNo.:

HK5.3

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

Unkonwn

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data