

For Research Use Only

Anti-Mouse CD197 (4B12)

Catalog Number: 65158-1-Ig



Basic Information

Catalog Number:

65158-1-Ig

Size:

100ug, 1 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

12775

Full Name:

chemokine (C-C motif) receptor 7

Purification Method:

Purified by protein-G affinity chromatography

CloneNo.:

4B12

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 5 ul Coralite® Plus 647-conjugated Anti-Mouse CD3 (CL647-65077, Clone: 17A2), and 1 ug Anti-Mouse CD197 (65158-1-Ig, Clone: 4B12) and FITC Anti-Rat IgG2a Antibody at dilution 1:200. Cells were not fixed.



1X10⁶ mouse splenocytes were surface stained with 5 ul Coralite® Plus 647-conjugated Anti-Mouse CD3 (CL647-65077, Clone: 17A2), and 1 ug Rat IgG2a Isotype Control and FITC Anti-Rat IgG2a Antibody at dilution 1:200. Cells were not fixed.