

Anti-Human LRP1 (5A6)

Catalog Number: 65233-1-Ig

Basic Information

Catalog Number: 65233-1-Ig	GenBank Accession Number: BC045107	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 4035	CloneNo.: 5A6
Source: Mouse	UNIPROT ID: Q07954	
Isotype: IgG2b	Full Name: low density lipoprotein-related protein 1 (alpha-2-macroglobulin receptor)	
	Calculated MW: 505 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information


LRP1 (Prolow-density lipoprotein receptor-related protein 1), also known as A2MR, APOER and CD91, is a type I transmembrane protein that belongs to the LDLR family. LRP1 is synthesized as a transmembrane glycosylated precursor of 600 kDa, and is subsequently cleaved to generate a large 515-kDa N-terminal extracellular subunit and an 85-kDa C-terminal transmembrane subunit, which are noncovalently associated with one another (PMID: 2112085; 8546712). LRP1 is a multifunctional receptor involved in the clearance of a large number of diverse ligands and modulating various cellular processes. LRP1 is implicated in conditions such as atherosclerosis, cancer and neurodegenerative diseases (PMID: 28381441).

Storage

Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

Selected Validation Data



1X10⁶ U-87 MG cells were surface stained with 1 ug Anti-Human LRP1 (65233-1-Ig, Clone:5A6) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (red), or 1 ug Control Antibody. Cells were not fixed.