

For Research Use Only

Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody, PBS Only

Catalog Number: 65571-1-PBS



Basic Information

Catalog Number:

65571-1-PBS

Concentration:

1mg, 2 mg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Applications

Tested Applications:

FC

Species Specificity:

n/a

Background Information

Storage

Storage:

Store at -80°C. Stable for one year after shipment.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.2 ug Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody (65571-1-MR, Clone:MOPC-173) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.2 ug Mouse IgG2a Isotype Control (MOPC-173) (65244-1-Ig, Clone: MOPC-173) (blue), or 0.2 ug Anti-Human CD8 (UCHT4) (65204-1-Ig, Clone: UCHT4)(green). Cells were not fixed. Lymphocytes were gated. This data was



1x10⁶ human PBMCs were surface stained with 0.2 ug Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody, PBS Only (65571-1-PBS, Clone: MOPC-173) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1) (red), or 0.2 ug Mouse IgG2a Isotype Control (MOPC-173) (65244-1-Ig, Clone: MOPC-173) (blue), or 0.2 ug Anti-Human CD8 (UCHT4) (65204-1-Ig, Clone: UCHT4) (green). Cells were not fixed. Lymphocytes were gated.