

For Research Use Only

Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody



Catalog Number: 65607-1-MR

Basic Information

Catalog Number:

65607-1-MR

Size:

100ug, 1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

24796

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

W3/13

Applications

Tested Applications:

FC

Species Specificity:

Rat

Background Information

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat splenocytes were surface stained with 0.25 ug Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed.



1x10⁶ rat splenocytes were surface stained with 0.25 ug Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) (red) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) (blue) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed.