

For Research Use Only

# Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody, PBS Only

Catalog Number: 65626-1-PBS



## Basic Information

Catalog Number:

65626-1-PBS

Concentration:

1mg, 2000 ug/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

Full Name:

CD1c molecule

Purification Method:

Protein A purification

CloneNo.:

L161

## Applications

Tested Applications:

FC

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human peripheral blood lymphocytes were surface stained with 0.25 ug Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody (65626-1-MR, Clone: L161) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) and PE-conjugated Goat Anti-Mouse IgG. Cells were co-stained with CoraLite® Plus 647 Anti-Human CD19. Cells were not fixed. This data was developed using the same antibody clone with 65626-1-PBS in a different storage buffer formulation.