

For Research Use Only

V5-tag Monoclonal antibody

Catalog Number: 66007-1-Ig **33 Publications**



Basic Information

Catalog Number:

66007-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

1 kDa

Purification Method:

Protein G purification

CloneNo.:

2B7A9

Applications

Tested Applications:

ELISA

Cited Applications:

CoIP, IP, RIP, WB

Species Specificity:

recombinant protein

Cited Species:

human

Background Information

The antibody is developed against a synthetic 14 aa peptide (GKPIPNNLLGLDST), the V5 tag, which presents in many vectors.

Notable Publications

Author	Pubmed ID	Journal	Application
Jinbo Fan	32967926	Life Sci Alliance	WB
Xing Sun	30237395	Cell Death Dis	
Ke Zheng	31777585	J Cancer	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

