

For Research Use Only

TST Monoclonal antibody

Catalog Number: 66018-2-Ig



Basic Information

Catalog Number: 66018-2-Ig	GenBank Accession Number: BC010148	Purification Method: Protein A purification
Size: 447 µg/ml	GeneID (NCBI): 7263	CloneNo.: 1D3D7
Source: Mouse	UNIPROT ID: Q16762	
Isotype: IgG2b	Full Name: thiosulfate sulfurtransferase (rhodanese)	
Immunogen Catalog Number: AG9235	Calculated MW: 297 aa, 33 kDa	
	Observed MW: 30-33 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data