

For Research Use Only

Bcl-XL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66020-2-PBS



Basic Information

Catalog Number:	66020-2-PBS	GenBank Accession Number:	BC019307	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	598	CloneNo.:	2B1A4
Source:	Mouse	UNIPROT ID:	Q07817		
Isotype:	IgG2a	Full Name:	BCL2-like 1		
Immunogen Catalog Number:	AG18037	Calculated MW:	26 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50613-1, Bcl-xL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66020-2-PBS. Detection antibody: 66020-1-PBS. Standard:Ag18037. Range: 1.563-100 ng/mL