

For Research Use Only

NOB1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66048-3-PBS



Basic Information

Catalog Number: 66048-3-PBS	GenBank Accession Number: BC000050	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 28987	CloneNo.: 1H1F10
Isotype: IgG2a	UNIPROT ID: Q9ULX3	
Immunogen Catalog Number: AG17907	Full Name: NIN1/RPN12 binding protein 1 homolog (S. cerevisiae)	
	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51874-1, NOB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66048-1-PBS. Detection antibody: 66048-3-PBS.



Cytometric bead array standard curve of MP51874-1, NOB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66048-1-PBS. Detection antibody: 66048-3-PBS. Standard:Ag17907. Range: 0.781-100 ng/mL.