

For Research Use Only

# NDUFA4L2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66050-5-PBS



## Basic Information

Catalog Number:

66050-5-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9233

GenBank Accession Number:

BC011910

GeneID (NCBI):

56901

UNIPROT ID:

Q9NRX3

Full Name:

NADH dehydrogenase (ubiquinone) 1  
alpha subcomplex, 4-like 2

Calculated MW:

87 aa, 10 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3C1C5

## Applications

Tested Applications:

IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP50529-3, NDUFA4L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66050-4-PBS. Detection antibody: 66050-5-PBS. Standard:Ag9233. Range: 0.195-100 ng/mL