

For Research Use Only

PA2G4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66055-2-PBS

Basic Information

Catalog Number:

66055-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7850

GenBank Accession Number:

BC001951

GeneID (NCBI):

5036

UNIPROT ID:

Q9UQ80

Full Name:

proliferation-associated 2G4, 38kDa

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

3G12F3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50029-1, PA2G4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66055-2-PBS. Detection antibody: HRP-conjugated 66055-1-PBS. Standard: Ag7850. Range: 1.25-20 ng/mL.