

For Research Use Only

Fetuin-A/AHSG Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66094-2-PBS



Basic Information

Catalog Number:

66094-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG10073

GenBank Accession Number:

BC048198

GeneID (NCBI):

197

UNIPROT ID:

P02765

Full Name:

alpha-2-HS-glycoprotein

Calculated MW:

367 aa, 39 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F6G10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51213-1, Fetuin-A/AHSG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66094-2-PBS. Detection antibody: 66094-3-PBS. Standard:Ag10073. Range: 0.098-100 ng/mL