

For Research Use Only

AK2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66127-2-PBS



Basic Information

Catalog Number:

66127-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG17911

GenBank Accession Number:

BC009405

GeneID (NCBI):

204

UNIPROT ID:

P54819

Full Name:

adenylate kinase 2

Calculated MW:

26 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

5D4E1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50407-1, AK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66127-1-PBS. Detection antibody: 66127-2-PBS. Standard:Ag17911. Range: 0.195-50 ng/mL