

For Research Use Only

MAP1LC3A,LC3 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66139-4-PBS



Basic Information

Catalog Number:

66139-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17959

GenBank Accession Number:

BC015810

GeneID (NCBI):

84557

UNIPROT ID:

Q9H492

Full Name:

microtubule-associated protein 1

light chain 3 alpha

Calculated MW:

121 aa, 14 kDa

Purification Method:

Protein G purification

CloneNo.:

1B5H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51371-1, MAP1LC3A,LC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66139-4-PBS. Detection antibody: 66139-5-PBS.



Cytometric bead array standard curve of MP51371-1, MAP1LC3A,LC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66139-4-PBS. Detection antibody: 66139-5-PBS. Standard:Ag17959. Range: 0.391-400 ng/mL.