

For Research Use Only

# IL-6 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66146-2-PBS

## Basic Information

|                                      |  |  |
|--------------------------------------|--|--|
| Catalog Number:<br>66146-2-PBS       | GenBank Accession Number:<br>BC015511            | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>3569                           | CloneNo.:<br>3B11C11                           |
| Source:<br>Mouse                     | ENSEMBL Gene ID:<br>ENSG00000136244              |  |
| Isotype:<br>IgG2b                    | UNIPROT ID:<br>P05231                            |  |
| Immunogen Catalog Number:<br>AG10646 | Full Name:<br>interleukin 6 (interferon, beta 2) |  |
|                                      | Calculated MW:<br>212 aa, 24 kDa                 |  |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50032-2, IL-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68659-3-PBS. Detection antibody: 66146-2-PBS. Standard:HZ-1019 . Range: 1.563-100 ng/mL.