

For Research Use Only

IL-1F6 / IL-36 α Monoclonal antibody



Catalog Number: 66147-1-Ig

Basic Information

Catalog Number:

66147-1-Ig

Size:

990 μ g/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG21487

GenBank Accession Number:

BC107043

GeneID (NCBI):

27179

UNIPROT ID:

Q9UHA7

Full Name:

interleukin 1 family, member 6 (epsilon)

Calculated MW:

158 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

1F2B8

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 66147-1-Ig (IL-1F6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.