

For Research Use Only

IL-1F6 / IL-36 α Monoclonal antibody



Catalog Number: 66147-1-Ig

Basic Information

Catalog Number: 66147-1-Ig	GenBank Accession Number: BC107043	Purification Method: Protein A purification
Size: 990 μ g/ml	GeneID (NCBI): 27179	CloneNo.: 1F2B8
Source: Mouse	UNIPROT ID: Q9UHA7	
Isotype: IgG2b	Full Name: interleukin 1 family, member 6 (epsilon)	
Immunogen Catalog Number: AG21487	Calculated MW: 158 aa, 18 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 66147-1-Ig (IL-1F6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.