

Ceruloplasmin Monoclonal antibody

Catalog Number: 66156-1-Ig **6 Publications**

Basic Information

Catalog Number:

66156-1-Ig

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15535

GenBank Accession Number:

BC146663

GeneID (NCBI):

1356

UNIPROT ID:

P00450

Full Name:

ceruloplasmin (ferroxidase)

Calculated MW:

1065 aa, 122 kDa

Observed MW:

140 kDa

Purification Method:

Protein G purification

CloneNo.:

3C8A5

Recommended Dilutions:

WB: 1:1000-1:5000

IHC: 1:20-1:200

IF/ICC: 1:400-1:1600

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : human placenta tissue, rat liver tissue, HSC-T6 cells, human kidney tissue, human blood tissue

IHC : human ovary tumor tissue, human placenta tissue, human liver tissue

IF/ICC : HeLa cells,

Background Information

Ceruloplasmin (CP) is a serum glycoprotein synthesized by hepatocytes and functions as the major transport protein of copper in blood where it acts as ferroxidase and is required for iron transport of transferrin. The intact ceruloplasmin migrated around 130-140 kDa and some proteolytic cleaved fragments could also be observed.

Notable Publications

Author	Pubmed ID	Journal	Application
Ni Yang	31706944	J Proteomics	WB
Marta Zarà	34360827	Int J Mol Sci	WB
Yong Zhang	34449964	FEBS Open Bio	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




rat liver tissue were subjected to SDS PAGE followed by western blot with 66156-1-Ig (Ceruloplasmin antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.




human placenta tissue were subjected to SDS PAGE followed by western blot with 66156-1-Ig (Ceruloplasmin antibody at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 66156-1-Ig (Ceruloplasmin antibody at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using 66156-1-Ig (Ceruloplasmin antibody at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using Ceruloplasmin antibody (66156-1-Ig, Clone: 3C8A5) at dilution of 1:800 and Coralite®488-Conjugated Goat Anti-Mouse IgG(H+L).