

For Research Use Only

BMI1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66161-2-PBS



Basic Information

Catalog Number:

66161-2-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG21284

GenBank Accession Number:

BC011652

GeneID (NCBI):

648

UNIPROT ID:

P35226

Full Name:

BMI1 polycomb ring finger oncogene

Calculated MW:

37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E3E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51443-1, BMI1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66161-2-PBS. Detection antibody: 66161-1-PBS.



Cytometric bead array standard curve of MP51443-1, BMI1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66161-2-PBS. Detection antibody: 66161-1-PBS. Standard: Ag21284. Range: 0.391-100 ng/mL.