

For Research Use Only



Vitamin D binding protein Monoclonal antibody, PBS Only

Catalog Number: 66175-1-PBS

Basic Information

Catalog Number: 66175-1-PBS	GenBank Accession Number: BC057228	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 2638	CloneNo.: 1E4D10
Source: Mouse	UNIPROT ID: P02774	
Isotype: IgG1	Full Name: group-specific component (vitamin D binding protein)	
Immunogen Catalog Number: AG9803	Calculated MW: 474 aa, 53 kDa	
	Observed MW: 52-58 kDa	

Applications

Tested Applications:
Indirect ELISA, IHC, WB

Species Specificity:
human

Background Information

Vitamin D binding protein is a sparsely glycosylated serum protein responsible for highly specific binding and tissue-specific delivery of vitamin D and its metabolites. In addition, it is also an actin scavenger, and is the precursor to the immunomodulatory protein, Gc-MAF. Vitamin D binding protein has been proposed to have significant roles in C5a chemotaxis, osteoclast development and possibly in macrophage activation/recruitment.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




human testis tissue were subjected to SDS PAGE followed by western blot with 66175-1-Ig (Vitamin D binding protein antibody at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66175-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 66175-1-Ig (Vitamin D binding protein antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66175-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 66175-1-Ig (Vitamin D binding protein antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66175-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 66175-1-Ig (Vitamin D binding protein antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 66175-1-PBS in a different storage buffer formulation.



1X10⁶ U-937 cells were intracellularly stained with 0.4 ug Anti-Human Vitamin D binding protein (66175-1-Ig, Clone:1E4D10) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 66175-1-PBS in a different storage buffer formulation.