

For Research Use Only

# P62,SQSTM1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66184-3-PBS



## Basic Information

Catalog Number:

66184-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13131

GenBank Accession Number:

BC017222

GeneID (NCBI):

8878

UNIPROT ID:

Q13501

Full Name:

sequestosome 1

Calculated MW:

48 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3E9C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51211-1, P62;SQSTM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66184-2-PBS. Detection antibody: 66184-3-PBS. Standard:Ag13131. Range: 0.781-100 ng/mL