

For Research Use Only

VAPB Monoclonal antibody

Catalog Number: 66191-1-Ig

Featured Product

22 Publications



Basic Information

Catalog Number:

66191-1-Ig

Concentration:

1300 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6313

GenBank Accession Number:

BC001712

GeneID (NCBI):

9217

UNIPROT ID:

O95292

Full Name:

VAMP (vesicle-associated membrane protein)-associated protein B and C

Calculated MW:

27 kDa

Observed MW:

30-32 kDa

Purification Method:

Protein A purification

CloneNo.:

4F6A6

Recommended Dilutions:

WB 1:1000-1:6000

IHC 1:20-1:200

IF/ICC 1:400-1:1600

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IF, CoIP

Species Specificity:

human

Cited Species:

human, mouse, monkey, yeast

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : A549 cells, HepG2 cells, LNCaP cells, MCF-7 cells, HeLa cells, Jurkat cells, K-562 cells

IHC : human breast cancer tissue, mouse kidney tissue

IF/ICC : HepG2 cells,

Background Information

Vesicle-associated membrane protein-associated protein B (VAPB) is an integral membrane protein localized to the endoplasmic reticulum (ER) membrane. VAPB has been implicated in various cellular processes, including ER stress, the unfolded protein response (UPR) and calcium homeostasis regulation. The mutations in the gene of VAPB cause amyotrophic lateral sclerosis 8 (ALS8) and some other related forms of motor neuron disease including late onset spinal muscular atrophy.

Notable Publications

Author	Pubmed ID	Journal	Application
Qiaoxia Zheng	36198318	Cell	WB
Christina James	31519755	J Biol Chem	WB
Alexa Bishop	31198905	Contact (Thousand Oaks)	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of VAPB antibody (66191-1-Ig, 1:5000) with si-control and si-VAPB transfected HepG2 cells.




Various lysates were subjected to SDS PAGE followed by western blot with 66191-1-Ig (VAPB antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using VAPB antibody (66191-1-Ig, Clone: 4F6A6) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 66191-1-Ig (VAPB Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 66191-1-Ig (VAPB Antibody) at dilution of 1:50 (under 40x lens).