

For Research Use Only

NINJ2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66202-4-PBS



Basic Information

Catalog Number:

66202-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5343

GenBank Accession Number:

BC057766

GeneID (NCBI):

4815

UNIPROT ID:

Q9NZG7

Full Name:

ninjurin 2

Calculated MW:

142 aa, 16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5E8C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50693-2, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-4-PBS. Detection antibody: 66202-3-PBS. Standard:Ag5343. Range: 0.391-100 ng/mL.