

For Research Use Only

# UCH-L1/PGP9.5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66230-3-PBS



## Basic Information

Catalog Number:

66230-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6547

GenBank Accession Number:

BC000332

GeneID (NCBI):

7345

UNIPROT ID:

P09936

Full Name:

ubiquitin carboxyl-terminal esterase

L1 (ubiquitin thiolesterase)

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D11H10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50875-1, UCHL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66230-3-PBS. Detection antibody: 66230-4-PBS. Standard:Ag6547. Range: 1.563-100 ng/mL.