

For Research Use Only

CD68 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66231-3-PBS



Basic Information

Catalog Number:

66231-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22815

GenBank Accession Number:

BC015557

GeneID (NCBI):

968

UNIPROT ID:

P34810

Full Name:

CD68 molecule

Calculated MW:

37 kDa

Purification Method:

Protein G purification

CloneNo.:

3E12D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50230-1, CD68 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66231-1-PBS. Detection antibody: 66231-3-PBS. Standard:Ag22815. Range: 0.391-25 ng/mL