

For Research Use Only

# PHOX2B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66254-4-PBS



## Basic Information

Catalog Number:

66254-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13196

GenBank Accession Number:

BC017199

GeneID (NCBI):

8929

UNIPROT ID:

Q99453

Full Name:

paired-like homeobox 2b

Calculated MW:

31 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A2B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51004-2, PHOX2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66254-4-PBS. Detection antibody: 66254-5-PBS. Standard:Ag13196. Range: 0.098-100 ng/mL.