

For Research Use Only

PHOX2B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66254-6-PBS



Basic Information

Catalog Number:

66254-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13196

GenBank Accession Number:

BC017199

GeneID (NCBI):

8929

UNIPROT ID:

Q99453

Full Name:

paired-like homeobox 2b

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

1F4B2

Applications

Tested Applications:

IHC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51004-3, PHOX2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66254-6-PBS. Detection antibody: 66254-3-PBS. Standard:Ag13196. Range: 0.195-100 ng/mL.