

For Research Use Only

OTUD7B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66276-2-PBS



Basic Information

Catalog Number:

66276-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10513

GenBank Accession Number:

BC020622

GeneID (NCBI):

56957

UNIPROT ID:

Q6GQQ9

Full Name:

OTU domain containing 7B

Calculated MW:

427aa, 48 kDa; 843aa, 93 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H6C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51475-1, OTUD7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66276-1-PBS. Detection antibody: 66276-2-PBS.



Cytometric bead array standard curve of MP51475-1, OTUD7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66276-1-PBS. Detection antibody: 66276-2-PBS. Standard:Ag10513. Range: 0.781-100 ng/mL