

For Research Use Only

SDCCAG8 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66284-3-PBS

Basic Information

Catalog Number: 66284-3-PBS	GenBank Accession Number: BC032454	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 10806	CloneNo.: 1F6G8
Source: Mouse	UNIPROT ID: Q86SQ7	
Isotype: IgG1	Full Name: serologically defined colon cancer antigen 8	
Immunogen Catalog Number: AG4264	Calculated MW: 713 aa, 83 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50102-1, SDCCAG8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66284-2-PBS. Detection antibody: 66284-3-PBS. Standard:Ag4264. Range: 0.391-100 ng/mL