

For Research Use Only

OGDH Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66285-3-PBS



Basic Information

Catalog Number:

66285-3-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7484

GenBank Accession Number:

BC004964

GeneID (NCBI):

4967

UNIPROT ID:

Q02218

Full Name:

oxoglutarate (alpha-ketoglutarate) dehydrogenase (lipoamide)

Calculated MW:

116 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G11E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51450-1, OGDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66285-2-PBS. Detection antibody: 66285-3-PBS.



Cytometric bead array standard curve of MP51450-1, OGDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66285-2-PBS. Detection antibody: 66285-3-PBS. Standard: Ag7484. Range: 0.098-100 ng/mL.