

For Research Use Only

# EPCAM/CD326 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66316-4-PBS



## Basic Information

Catalog Number:

66316-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15393

GenBank Accession Number:

BC014785

GeneID (NCBI):

4072

UNIPROT ID:

P16422

Full Name:

epithelial cell adhesion molecule

Calculated MW:

314 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

3D2E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50309-2, EPCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66316-4-PBS. Detection antibody: 66316-3-PBS. Standard:Ag15393. Range: 0.781-25 ng/mL