

For Research Use Only

EPCAM/CD326 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66316-6-PBS



Basic Information

Catalog Number:

66316-6-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG15393

GenBank Accession Number:

BC014785

GeneID (NCBI):

4072

UNIPROT ID:

P16422

Full Name:

epithelial cell adhesion molecule

Calculated MW:

314 aa, 35 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3G10C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50309-4, EPCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66316-6-PBS. Detection antibody: 66316-7-PBS. Standard:Ag15393. Range: 3.125-100 ng/mL.