

For Research Use Only

BCL6 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66340-3-PBS

Basic Information

Catalog Number:

66340-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15519

GenBank Accession Number:

BC150184

GeneID (NCBI):

604

UNIPROT ID:

P41182

Full Name:

B-cell CLL/lymphoma 6

Calculated MW:

706 aa, 79 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H4D2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80008-2, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-3-PBS. Detection antibody: 81756-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL