

For Research Use Only

MKP-2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66349-2-PBS

Basic Information

Catalog Number: 66349-2-PBS	GenBank Accession Number: BC014565	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 1846	CloneNo.: 1G12G6
Source: Mouse	UNIPROT ID: Q13115	
Isotype: IgG2a	Full Name: dual specificity phosphatase 4	
Immunogen Catalog Number: AG1206	Calculated MW: 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50101-1, MKP-2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66349-2-PBS. Detection antibody: 66349-3-PBS. Standard:Ag1206. Range:6.25-100 ng/mL.