

For Research Use Only

# MKP-2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66349-3-PBS

## Basic Information

Catalog Number: 66349-3-PBS	GenBank Accession Number: BC014565	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 1846	CloneNo.: 1D10C12
Source: Mouse	UNIPROT ID: Q13115	
Isotype: IgG1	Full Name: dual specificity phosphatase 4	
Immunogen Catalog Number: AG1206	Calculated MW: 43 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50101-1, MKP-2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66349-2-PBS. Detection antibody: 66349-3-PBS. Standard:Ag1206. Range:6.25-100 ng/mL.