

For Research Use Only

# EBI3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66354-2-PBS



## Basic Information

Catalog Number:

66354-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3051

GenBank Accession Number:

BC015364

GeneID (NCBI):

10148

UNIPROT ID:

Q14213

Full Name:

Epstein-Barr virus induced 3

Calculated MW:

229 aa, 25 kDa

Purification Method:

Protein G purification

CloneNo.:

1A1G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51311-1, EBI3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66354-2-PBS. Detection antibody: 66354-3-PBS. Standard:Ag3051. Range: 0.098-100 ng/mL.