

For Research Use Only

Mouse IgG3 isotype control Monoclonal antibody



Catalog Number: 66360-4-Ig

1 Publications

Basic Information

Catalog Number:

66360-4-Ig

Source:

Mouse

Isotype:

IgG3

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

1H4A5

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

RIP

Species Specificity:

megathura crenulata

Positive Controls:

WB :

Background Information

This antibody is used as the isotype control of mouse IgG3. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Notable Publications

Author	Pubmed ID	Journal	Application
Youjun Li	40294852	Biochim Biophys Acta Mol Basis Dis	RIP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66360-4-Ig (KLH antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.