

For Research Use Only

MMP2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66366-3-PBS



Basic Information

Catalog Number:

66366-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25039

GenBank Accession Number:

BC002576

GeneID (NCBI):

4313

UNIPROT ID:

P08253

Full Name:

matrix metalloproteinase 2
(gelatinase A, 72kDa gelatinase,
72kDa type IV collagenase)

Calculated MW:

72 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H10G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51054-1, MMP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66366-2-PBS. Detection antibody: 66366-3-PBS. Standard:Ag25039. Range: 0.098-100 ng/mL.