

For Research Use Only

NFIL3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66368-2-PBS



Basic Information

Catalog Number: 66368-2-PBS	GenBank Accession Number: BC008197	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 4783	CloneNo.: 2D6A6
Isotype: IgG1	UNIPROT ID: Q16649	
Immunogen Catalog Number: AG2380	Full Name: nuclear factor, interleukin 3 regulated	
	Calculated MW: 462 aa, 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51436-1, NFIL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66368-2-PBS. Detection antibody: 66368-3-PBS.



Cytometric bead array standard curve of MP51436-1, NFIL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66368-2-PBS. Detection antibody: 66368-3-PBS. Standard: Ag2380. Range: 3.125-400 ng/mL.