

For Research Use Only

PLA2G1B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66397-2-PBS



Basic Information

Catalog Number:

66397-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8812

GenBank Accession Number:

BC005386

GeneID (NCBI):

5319

UNIPROT ID:

P04054

Full Name:

phospholipase A2, group IB (pancreas)

Calculated MW:

70aa, 8 kDa; 148aa, 16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E6D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51376-1, PLA2G1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66397-1-PBS. Detection antibody: 66397-2-PBS.



Cytometric bead array standard curve of MP51376-1, PLA2G1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66397-1-PBS. Detection antibody: 66397-2-PBS. Standard:Ag8812. Range: 0.391-400 ng/mL