

For Research Use Only

NFKBIA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66418-3-PBS



Basic Information

Catalog Number:

66418-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22040

GenBank Accession Number:

BC004983

GeneID (NCBI):

4792

UNIPROT ID:

P25963

Full Name:

nuclear factor of kappa light

polypeptide gene enhancer in B-cells

inhibitor, alpha

Calculated MW:

36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B11D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51027-1, NFKBIA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66418-2-PBS. Detection antibody: 66418-3-PBS. Standard:Ag0345. Range: 0.391-50 ng/mL