

For Research Use Only

CDC37 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66420-3-PBS



Basic Information

Catalog Number: 66420-3-PBS	GenBank Accession Number: BC000083	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 11140	CloneNo.: 4B5B12
Isotype: IgG1	UNIPROT ID: Q16543	
Immunogen Catalog Number: AG0267	Full Name: cell division cycle 37 homolog (S. cerevisiae)	
	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51702-1, CDC37 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66420-3-PBS. Detection antibody: 66420-2-PBS.



Cytometric bead array standard curve of MP51702-1, CDC37 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66420-3-PBS. Detection antibody: 66420-2-PBS. Standard:Ag0267. Range: 1.563-200 ng/mL.